



## KODAK Direct View CR Elite System

System

KODAK DirectView Elite CR

imaging),

(long-length

# KODAK Direct View CR Elite System

▶ 35 84 / 15 30

▶ (body parts)

▶ KODAK DirectView Elite CR

▶ On-line IT- - CSA (Cisco Security Agent)

▶ DICOM 3.0 (DICOM Print SCU, DICOM Storage SCU, DICOM Worklist SCU)

▶ Kodak DirectView CR HR, EHR-M2, portable vertical long-length, PQ


▶ 122

▶ 33

▶ 90

▶ 79

▶ 10 / 20 / 16 / 12 /



KODAK DirectView CR Carestream Health.

KODAK DirectView  
ROP(Remote Operation Panel)

KODAK DirectView Captur Link

KODAK DirectView CR

CR (DR) KODAK  
DirectView

EHR-M2,

CR

5 ROP  
CR KODAK DirectView.

ROP.

EHR-M2

ROP 8

KODAK DirectView TQT

48

R.

TQT

10

KODAK DirectView EVP Plus

KODAK DirectView CR

LLI.

43 129

Kodak Long Length  
Imaging.

▶EVP Plus

▶EVP Plus

35 84 )

▶EVP Plus



# KODAK Direct View CR Elite System

10	/	
20	/	
- 33		
18	x 24	100
24	x 30	76
14 in.	x 14 in./35 x 35	102*
14 in.	x 17in./35 x 43	90*
15	x 30 (Dental)	122
35	x 84 (Long Length)	88
*		

Kodak DirectView CR cassette / with PQ Rigid screen / 35 x 43 c , 35 x 84 c , 35 x 35 , 24 x 30 , 18 x 24 , 15 x 30

Kodak DirectView CR cassette / with PQ Rigid screen / 35 x 43 35 x 35 ( 10 / )

Kodak DirectView CR cassette / with HR screen / 24 x 30 18 x 24 ( 10 / )

Kodak DirectView CR mammography cassette / with EHR-M2 screen / 24 x 30 18 x 24 ( 20 / )

65536 (16 / )

: 50

4096 (12 / )

: 65

: 105

: 140

10 BaseT / 100 BaseT

12

220-240 V, 50/60 , 5

CarestreamHealth

199991 , ,6, -140

(495) 236 44 57  
(495) 236 9814

(495)545 3538

E-mail:  
Info@neorentgen.com

<http://www.neorentgen.com>